



Teeth Whitening Strips

User manual

Warning:

1. Not intended for children under 16 years of age.
2. Do not swallow the strip.
3. Avoid contact of the strips with eyes or skin.
4. If the strip comes into contact with the skin, it may cause temporary skin irritation.
5. Stop using and consult a dentist immediately if any discomfort occurs.
6. Not suitable for people with ulcers, exposed gums, or after oral surgery.

Instructions for Use:

1. Remove the strips from the backing.
2. Apply the upper strip to the upper teeth, and repeat this with the lower teeth.
3. Leave on for 30 minutes; avoid eating or drinking.
4. Remove and dispose of the strips, rinse off any remaining gel if necessary, and admire the result.

Shelf Life:

2 years

Storage:

Store in a dry and cool place (keep out of reach of young children).

Disposal:

Dispose of the product after use in accordance with local regulations.

Important Tips:

1. Use the teeth whitening strips regularly, every three months, for a whiter and brighter smile.
2. Reduce or avoid smoking, and drinking coffee, tea, cola, or red wine, as well as eating strongly colored foods such as cherries and blueberries to achieve a better and longer-lasting effect.
3. Use all strips within two weeks for the best result.

User Guide:

1. Avoid using the strips immediately after brushing your teeth. Wait at least one hour after brushing before using the strips.
2. Do not reuse the strips; dispose of them after one use.
3. The ingredients of the strips are safe; if a small amount of gel is swallowed, it is completely harmless.
4. Remove any remaining gel by brushing and rinsing your teeth.