



Teeth whitening kit

HOW DOES IT WORK?

Our Easy Smile whitener uses innovative blue LED light technology. The blue light reacts with the whitening gel to effectively remove stains from the surface of the teeth. The mouthpiece has been engineered to be as soft and comfortable as possible.

The Easy Smile Kit is comfortable to use and causes no discomfort during the whitening process. When you bite down on the mouthpiece, it fits snugly and fully covers your teeth. Before beginning the treatment, it is important to understand how it works and what steps to follow. Please read the instructions carefully before use.

APPLICATIONS

The whitener is ideal for removing:

Coffee stains, tobacco stains, genetic discoloration, tetracycline stains, food dyes, aging stains, and other causes of discoloration.

STORAGE

Do not expose to direct sunlight. Store in a cool place.

USER GUIDE

1. Rinse the mouthpiece and your mouth thoroughly before use.
2. Look in the mirror and compare your teeth to the color guide to determine your starting shade.
3. Twist the whitening gel pen until gel appears in the brush. Apply a thin, even layer of gel to all your teeth—both upper and lower.
Do not use too much gel (recommended amount is 0.5 ml – 0.7 ml per use). More gel will not improve or speed up results.
Avoid getting gel on the gums or inside of the lips. If applied accidentally, do not worry—any burning sensation will disappear within a few minutes.

Rinse with water afterwards.

4. Ensure the device is fully charged before the first use.
Place the mouthpiece in your mouth, bite down, and hold the button to activate the LED light. Breathe deeply to reduce saliva production (excess saliva dilutes the gel).
5. The device will automatically switch off after 16 minutes. Remove the mouthpiece, rinse both it and your mouth with water, and compare your teeth to the guide again.
6. Do not eat or drink anything except water for 30 minutes after treatment.

VISUAL INSTRUCTIONS (illustrations):



Rinse the mouth
before use



Apply the gel with the pen
on all teeth



Insert the mouthpiece and bite
down on it, then switch on the LED
light



Wait 16 minutes to
conclude the treatment



Rinse the mouth
with water



Wash mouthpiece with water
and close gel pen with cap

REMINDERS

1. Rinse your mouth thoroughly, dry the teeth completely, and apply the gel evenly to the surface of the teeth.
2. Excess saliva may dilute the effectiveness of the gel.
3. Everyone's teeth are different; some people may experience slight sensitivity — this is normal. A desensitizing gel can be applied to relieve discomfort. If sensitivity is severe, reduce usage time or discontinue use immediately.
4. Avoid applying the whitening gel directly to the gums or the inside of the lips.
5. If the gel accidentally gets on your skin or clothing, rinse it off immediately.
6. Always close the pen with the cap after each use to prevent the active ingredient from evaporating.
7. Store the product out of reach of children.
8. It is recommended not to use the whitener more than 4 times a week.

INGREDIENTS

Sorbitol, Water, Phthalimidoperoxycaproic acid (PAP), Propylene Glycol, Glycerin, Potassium Nitrate, PEG-8, Hydroxyapatite, Cellulose Gum, Hydroxyethylcellulose, Xanthan Gum, Sodium Saccharin, Menthol, Methylparaben, Aloe Barbadensis (Aloe Vera) Leaf Extract, Chamomile Extract, Pomegranate Extract, Sodium Bicarbonate, Propylparaben.

NOTES

- Charging time is approximately 2 hours.
- Charging status: the Easy Smile mouthpiece will flash while charging.
- Press and hold the button for 2 seconds to turn the device on or off.
- The expected lifespan of the whitener is approximately 24 months.

